

Introducing...

BRASSELER^{VE}TM

VISION ENHANCEMENT

*Innovative and Ergonomic Vision Enhancement Loupes
from the Leader in Dental Instrumentation*

Superior Resolution

Highest quality optical glass ensures perfect convergence and maximum image resolution, reducing eye fatigue.

Enhanced Field of View

Our large fields of view, with edge-to-edge clarity, ease your adjustment to utilizing enhanced vision in your practice.

Outstanding Depth of Field

Our quality optics ensures the best possible depth of field at all working distances. BrasselerVE loupes give you greater freedom, comfort and efficiency.

Professional Comfort & Style

BrasselerVE loupes offer a wide selection of frame types and colors in both lightweight titanium and popular sport styles.



Convenient & Personalized Service

Brasseler USA representatives are available nationwide to give you the in-office personal attention and service you require. Need your personal vision prescription in your loupes? Let our highly trained staff assist you.

Extensive Selection

BrasselerVE options include a wide range of magnifications and working distances in both flip-up and through-the-lens styles. For your advanced magnification needs, prism loupes are available in 4.0x and 5.0x flip-up styles.

Ergonomic Design

BrasselerVE loupes are ideally balanced and ultra-light, reducing fatigue. Plus, our flip-up style incorporates a unique double hinge system which allows you to customize them to your best viewing angle.

Risk Free Purchase

If you are not completely satisfied with your loupes or headlight purchase, simply return it within 30 days for a full refund.*

*Restrictions apply for prescription-added loupes.



Expect the Best. Buy Direct.

To order call 800.841.4522 or fax 888.610.1937.
Visit our web site: www.BrasselerUSA.com.



BRASSELER USA
DENTAL INSTRUMENTATION

BRASSELER **VE**TM

VISION ENHANCEMENT

Portable LED Headlight



Traditional LED Light



Fiber Optic Light



BrasselerVETM Portable LED Headlight

- *Battery pack features charge indicator lights, so you'll always know when it's time for the battery to be recharged*
- *Light intensity does not degrade as the battery charge decreases*

- *Powerful 5,000+ foot candle illumination*
- *Unique lens design ensures light beam remains focused, minimizing diffusion and scattering*
- *Pure white light with fully adjustable beam intensity prevents color distortion and glare*
- *Your BrasselerVE light can be customized to mount onto most competitive loupes or can be worn on our convenient headband*
- *Weighing less than one ounce, the feather light BrasselerVE headlight delivers supreme comfort and freedom*
- *State-of-the-art lithium polymer battery provides over 6 hours use-life per charge*



- *Optional composite filter available to prevent curing of composite materials*

Expect the Best. Buy Direct.

To order call **800.841.4522**, in Canada **800.363.3838**.



BRASSELER USA[®]
DENTAL INSTRUMENTATION